

# IGEL

## BASE CURVE SELECTION GUIDE

<b>DIAMETER (mm)</b>	<b>13.50</b>	<b>13.80</b>	<b>14.00</b>	<b>14.20</b>	<b>14.50</b>
<b>MEAN K (mm)</b>	<b>BASE CURVE (mm)</b>				
<b>7.20</b>			8.00	8.10	8.20
<b>7.30</b>		8.00	8.10	8.20	8.30
<b>7.40</b>	8.00	8.10	8.20	8.30	8.40
<b>7.50</b>	8.10	8.20	8.30	8.40	8.50
<b>7.60</b>	8.20	8.30	8.40	8.50	8.60
<b>7.70</b>	8.30	8.40	8.50	8.60	8.70
<b>7.80</b>	8.40	8.50	8.60	8.70	8.80
<b>7.90</b>	8.50	8.60	8.70	8.80	8.90
<b>8.00</b>	8.60	8.70	8.80	8.90	9.00
<b>8.10</b>	8.70	8.80	8.90	9.00	9.10
<b>8.20</b>	8.80	8.90	9.00	9.10	9.20
<b>8.30</b>	8.90	9.00	9.10	9.20	
<b>8.40</b>	9.00	9.10	9.20		