

IGEL

DIOPTERS TO RADIUS CONVERSION CHART

DIOPTERS D	RADIUS mm
39.00	8.65
39.25	8.60
39.50	8.54
39.75	8.49
40.00	8.44
40.25	8.39
40.50	8.33
40.75	8.28
41.00	8.23
41.25	8.18
41.50	8.13
41.75	8.08
42.00	8.04
42.25	7.99
42.50	7.94
42.75	7.89
43.00	7.85
43.25	7.80
43.50	7.76
43.75	7.71
44.00	7.67
44.25	7.63
44.50	7.58
44.75	7.54
45.00	7.50
45.25	7.46
45.50	7.42
45.75	7.38
46.00	7.34
46.25	7.30
46.50	7.26
46.75	7.22
47.00	7.18
47.25	7.14
47.50	7.11
47.75	7.07
48.00	7.03
48.25	6.99
48.50	6.96
48.75	6.92
49.00	6.89
49.25	6.85
49.50	6.82
49.75	6.78
50.00	6.75